

SPECTRA OPTIA®

APHERESIS SYSTEM

The Confidence You Expect with
Spectra Optia Apheresis System Enhancements

System Comparison

COBE® Spectra Apheresis System	Spectra Optia System
Manual Interface Management	Automated Interface Management with Operator Control
Mobile	Highly Mobile and Transportable
Manual Data Capture	Prints and Exports Data
Manual Data Retrieval	Access Historical Data
Push Button Panel, LED Display	Touch Screen, Graphical User Interface Display
Multiple Screen Navigation Required to Access Desired Data	Minimal Screen Navigation Required to Access Desired Data
One Tubing Set Type = One Procedure Type	One Tubing Set Type = Multiple Procedures Types
Multiple Steps Required to Load Tubing Set	Minimal Steps Required to Load Tubing Set
Involved Operator Training	Simplified Operator Training
Noticeable Operating Noise	Low Operating Noise
Stand Alone Tube Sealer	Incorporated Tube Sealer

Your next generation therapeutic apheresis and cell therapy platform delivers COBE Spectra system confidence – designed to enhance versatility, performance, and patient comfort and safety.

Terumo BCT, Inc.

10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA

USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT Europe N.V.

Europe, Middle East and Africa
Ikaroslaan 41
1930 Zaventem
Belgium

Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

Terumo BCT (Asia Pacific) Ltd.

Room 3903-3903A, 39/F
ACE Tower, Windsor House
311 Gloucester Road
Causeway Bay, Hong Kong

Phone: +852.2283.0700
Fax: +852.2576.1311

Terumo BCT Latin America S.A.

La Pampa 1517 - 12th Floor
C1428D2E
Buenos Aires
Argentina

Phone: + 54.11.5530.5200
Fax: + 54.11.5530.5201

Terumo BCT Japan, Inc.

Ebisu South One Bldg. 9F
1-7-8, Ebisu-minami,
Shibuya-ku, Tokyo 150-0022
Japan

Phone: +81.3.6743.7890
Fax: +81.3.6743.9800